

 BC Ambulance Service BC Emergency Health Services	110 VOLT TRANSFER SWITCH INSTALATION	November 6, 2017 Page 1 of 8
--	---	---------------------------------

AMBULANCE TYPE: CRESTLINE SINGLE STRETCHER

UNITS AFFECTED: ALL CRESTLINE SINGLE STRETCHER UNITS EQUIPPED
WITH THE STRYKER POWERPRO STRETCHER

ISSUE: To provide 110V power to the internal receptacles of the ambulance from shoreline power or the Go Power inverter, to facilitate charging of the Stryker PowerPro Stretcher battery.

ACTION: INSTALLATION OF AN 110V TRANSFER SWITCH

LABOUR TIME: 1.5 hrs **LABOUR CODE:** GMOTS

SERVICE PROCEDURE:

Please read through this modification entirely before starting modification.

Contact BCAS Fleet Operations if there are any questions or concerns.

1-877-652-7465

Check contents of kit provided.

- 2 Marrets - 110V connectors
- 1 female 110V plug
- 6 Self tapping screws #10 x 5/8
- 2 Cover Screws #10-32 x 5/8
- 2 cable clamps
- 3 Zap straps
- Mounting Bracket – Crestline
- Transfer Switch
- Plastic hole plug



BEFORE STARTING THIS MODIFICATION:

Confirm 110 volt power is not connected to the ambulance shoreline power and the engine is switched off.

110 VOLT TRANSFER SWITCH INSTALATION

November 6, 2017
Page 2 of 8


LOCATE THE 110 VOLT BREAKER BOX IN THE ELECTRICAL COMPARTMENT

STEP 1:

Remove the cover from the 110V breaker box and compare the wiring with the diagram below.

IMPORTANT!

Confirm existing wiring configuration before continuing.






110 VOLT TRANSFER SWITCH INSTALATION

STEP 2

1. Temporarily install the supplied mounting plate in the location shown. Use the breaker box faceplate cover screws to secure the mounting plate.



2. Using the mounting plate as a template, drill two 1/8" pilot holes in the wall and install the supplied mounting screws.

3. Install the Transfer Switch to the mounting plate using the supplied four screws.






110 VOLT TRANSFER SWITCH INSTALATION

STEP 3

1. Remove the knock out from the right center of the breaker box & install the cable clamp.
2. Measures 5 inches from the end of the transfer switch harness labeled "Connect to AC Load Only" and remove the outer harness insulation.
3. Feed the transfer switch harness into the clamp allowing up to ½ inch of outer harness insulation through the clamp and secure.



4. Measure 8 inches from the end of the transfer switch harness labeled "Connect to AC power source" and remove the outer case insulation.
5. Remove the knock out from the bottom left of the box & install the supplied cable clamp.
6. Feed the transfer switch harness into the clamp allowing ½ inch of outer harness insulation through the clamp and secure.



November 6, 2017
Page 5 of 8

Note: The dashed lines show the original wire path. The faded lines show the wiring not disturbed by this modification.

CHARGER

1. Remove this ground bus bar. Connect the green ground wires from the transfer switch harnesses "AC load & AC power" to the grounding bus bar, then re-attach the bus bar to the breaker box.

2. Remove the black wire from the 110v outlet breaker and connect using a Marret to the black wire from the AC load transfer switch.

3. Remove the 110v outlet white wire from the neutral bus bar and connect with a Marret to the white wire from the transfer switch "AC load" harness.

4. Connect the black wire from the transfer switch "AC power source" harness to the AC breaker.

5. Plug into Inverter power cord. See next page.

CONNECT TO AC POWER

FROM UPPER SHORE LINE PLUG

FROM LOWER SHORE LINE PLUG

CONNECT TO AC LOAD

TO A/C HEATER


TO 110V COMP OUTLETS

PROwatt SW TRANSFER SWITCH

CONNECT TO AC LOAD



110 VOLT TRANSFER SWITCH INSTALATION




STEP 4

Locate the power cord coming from the Go Power Inverter and going to the GFI plug in the patient compartment. Cut this cord being sure to leave enough length to allow for rerouting to connect to the "Transfer Switch Inverter Only" harness above the breaker box. Install the supplied female extension cord receptacle onto the Inverter cord and plug into the "Transfer Switch Inverter only" plug. Secure with supplied zap strap.




110 VOLT TRANSFER SWITCH INSTALATION

Step 5: Recheck all the new connections to match this diagram. Ensure all Marret connections are taped with electricians tape, and all screw connections are tight.



Step 6

Re-install the Breaker Box Cover and tighten all mounting screws.




Step 7

In the patient compartment, remove the lower box labeled "Inverter Power". Pull the power cord wire out with the box and recycle both. Install supplied hole plug with silicone.



Use a razor blade to trim and peel the remaining box decal to read:
Power 110 Volt AC



System Test:

Using a GFI protected circuit; plug a shoreline power cord into the lower 110V shore power receptacle on the outside of the ambulance.

Plug an 110V tester into each of the patient compartment 110 volt outlets – power should be present.
Note: There is a second outlet located in the kit tree.

Start the engine and turn the connect switch on. Unplug the shore power from the outside of the ambulance.

Confirm the 110V outlets in the patient compartment still have power.

Labour Code: GMOTS – use this labour code on the invoice and charge 1.5 hrs.